
RISE

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EVOS™ M5000 Imaging System (MICROSCOPIO A FLUORESCENZA)

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Valenzano (BA), 70010

ENG: Within the framework of the National Recovery and Resilience Plan (PNRR), the Department of Veterinary Medicine (DiMeV) of the University of Bari Aldo Moro contributes to the INF-ACT extended partnership (PE00000007), dedicated to emerging infectious diseases through an integrated One Health approach. DiMeV contributes particularly to Research Node 2, 'Arthropod Vectors and Vector-Borne Pathogens', which aims to deepen knowledge of vector bionomics and vector-host-pathogen interactions and to develop innovative approaches for surveillance and eco-friendly vector control.

What is this technique?

The EVOS M5000 Imaging System is a versatile digital inverted benchtop microscopy platform designed for a wide range of microscopy applications, helping to optimize the entire workflow from image acquisition to data analysis.

The main acquisition features include autofocus, Z-stack acquisition, time-lapse imaging, and multi-channel image acquisition. The integrated software allows users to directly manage

image acquisition, analysis, and processing through a range of tools, including brightness and contrast adjustment.

How it works?

The EVOS M5000 integrates precision optics, an 18.5-inch high-resolution articulated LCD monitor, and a 3.2 MP monochrome CMOS camera. The monochrome sensor provides high sensitivity for the detection of fluorescence signals, including low-intensity signals, making it particularly suitable for fluorescence imaging applications.

The system uses RGB LED illumination, enabling the acquisition of high-quality, low-noise color images using the same monochrome sensor. This configuration provides flexibility for both fluorescence and transmitted-light imaging applications.

The optical configuration supports imaging in three fluorescence channels, as well as visible-light imaging with phase contrast in transmitted light. The available objectives are 40X and 100X, providing high-magnification imaging capabilities for detailed sample observation.

Interchangeable sample holders allow the system to accommodate different sample types and vessels, including multi-well plates and cell-culture flasks, making it suitable for a variety of biological and cell-based applications.

The system requires no warm-up or cool-down periods, allowing image acquisition to begin rapidly and providing a streamlined and efficient imaging workflow.

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