RISE

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MALDI BIOTYPER SIRIUS

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Valenzano (BA), 70010

ENG: Within the framework of the PNRR, the University of Bari, in partnership with the University of Pavia and key stakeholders, participates in the INF-ACT research program, which was developed to address unmet needs related to emerging infectious diseases in humans. The program adopts a One Health approach, considering human health within a broader context that includes domestic and wild animals as potential disease reservoirs, as well as environmental factors that may increase the risk of spillover. The project is organized into five research nodes. Three “vertical” research nodes focus on specific research themes: **Node 1**, emerging and re-emerging viral diseases, with a focus on respiratory and zoonotic viruses; **Node 2**, arthropod vectors and vector-borne pathogens, with a focus on vector-borne diseases most at risk of expanding or emerging in Italy, such as arboviruses; and **Node 3**, diseases caused by bacteria and fungi resistant to multiple antibiotics, with a focus on antimicrobial resistance and the molecular mechanisms of multidrug resistance. In this context, the MALDI-TOF MS system — MALDI Biotyper® sirius System, Bruker Daltonics, Germany — was acquired in 2025 and assigned to the Bacteriology and Mycology

Laboratories of the Department of Veterinary Medicine. Thanks to this acquisition, the identification of bacteria and fungi from domestic and wild animals, as well as from environmental samples, has been simplified and improved. The purchase of the instrument has supported research conducted by PhD students, technical staff, faculty members, and students from Node 3 and other projects.

What is this technique?

MALDI-TOF MS stands for **Matrix-Assisted Laser Desorption/Ionization–Time of Flight Mass Spectrometry**. MALDI-TOF mass spectrometry is a technique used to rapidly identify microorganisms by analyzing their protein profiles. It generates characteristic spectra that act like molecular fingerprints for bacteria and fungi. In clinical microbiology, it is now widely used because it can identify microorganisms within minutes, replacing many slower culture-based identification methods.

How it works?

MALDI-TOF MS works by analyzing the protein profile of a microorganism. A fragment of biological material, such as part of a bacterial colony, is placed on a metal target plate and covered with a matrix solution. When a laser hits the sample, the matrix absorbs the energy, helping the microbial proteins vaporize and become electrically charged without being destroyed. These charged particles, or ions, are then accelerated through a vacuum tube called a flight tube. Smaller ions travel faster than larger ones, so they reach the detector at different times. The instrument records these times and produces a mass spectrum, which acts like a unique protein fingerprint. The software compares the fingerprint with reference spectra in a database containing thousands of profiles and quickly identifies the microorganism. In less than a minute, the algorithm calculates a score and identifies the pathogen at the genus and species level.

MALDI-TOF MS enhances research, diagnostics, and laboratory services by enabling the rapid and accurate identification of microorganisms, including those that are difficult to culture or phenotypically ambiguous, such as anaerobes, fungi, and non-tuberculous mycobacteria. Although the initial cost of the instrument is high, the cost per test is very low, making it more economical than many traditional biochemical identification methods. In our laboratory, MALDI-TOF MS has reduced bacterial identification times from 48 hours to just a few minutes, with accuracy reaching 99.9% under validated conditions, allowing faster selection of targeted antibiotic therapy and lowering operating costs.

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- Bruker — MALDI Biotyper® for Veterinary Microbiology. MALDI Biotyper® system as a rapid method for identifying cultivated veterinary microorganisms, including bacteria and fungi, supporting veterinary microbiology and zoonosis surveillance, from: <https://www.bruker.com/en/applications/microbiology-and-diagnostics/veterinary-microbiology/maldi-biotyper-for-veterinary.html>.