RISE

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Via Orabona 4, 70125 (BA)

ENG: The Orbitrap Exploris 480 is a high-resolution mass spectrometer designed for advanced proteomic and peptidomic applications, providing high mass accuracy, resolving power, sensitivity, and analytical robustness. In our configuration, the mass spectrometer is coupled to a Vanquish Neo UHPLC system and equipped with a nano-electrospray ionization (nano-ESI) source, enabling highly sensitive nanoLC-ESI-HRMS/MS analysis of complex peptide mixtures. In bottom-up proteomics workflows, proteins are enzymatically digested, typically using trypsin, and the resulting peptides are separated by reversed-phase nano-liquid chromatography prior to high-resolution tandem mass spectrometry analysis. Data-Dependent Acquisition (DDA) workflows can be applied, in which a high-resolution full MS scan is followed by the automatic selection of the most abundant precursor ions for fragmentation, typically by Higher-energy Collisional Dissociation (HCD), and subsequent MS/MS

analysis. This approach enables accurate peptide identification and comprehensive protein characterization. The system offers broad applications across food, microbiological, and environmental sciences. In food safety, proteomic approaches can support the characterization of foodborne pathogens, microbial stress responses, virulence-associated proteins, antimicrobial resistance mechanisms, allergens, and microbial interactions, contributing to the identification of potential molecular biomarkers of biological hazards. In environmental microbiology, differential proteomics can provide insights into microbial responses to environmental stressors, contaminants, and changing physicochemical conditions, supporting the assessment of microbial adaptation and ecosystem functionality. Peptidomic applications further enable the identification of bioactive peptides, (e.g., antimicrobial, antioxidant, and enzyme-inhibitory molecules) generated during fermentation, food processing, or gastrointestinal digestion.

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