
RISE

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Via Orabona 4, 70125 (BA)

ENG: Within the framework the PNRR (National Recovery and Resilience Plan) Agritech and MAR.V.E.L., the University of Bari has equipped DiSSPA with an integrated Gas Chromatography-Mass Spectrometry-Olfactometry system (GC-MS-ODP), combining an Agilent Gas Chromatograph-Mass Spectrometer with a PAL automated sampling system and an Olfactory Detection Port for the simultaneous instrumental and human-nose detection of volatile and odour-active compounds. The GC column effluent is split between the mass spectrometer, which enables identification and semi-quantification of volatile organic compounds (VOCs), and the heated ODP transfer line, at which trained assessors perceive, in real time, the intensity, duration and identify of the odour associated with each eluting

compound (if aroma-active), linking chemical data directly to its sensory perception. This methodology is particularly suited to the study of any type of food and ingredient. Combining GC-MS-ODP data with sensory panel evaluation and compositional analysis supports an integrated approach linking chemical composition, volatilome and perceived sensory quality. The objective is to provide this VOC and olfactometric characterisation service to public and private stakeholders requiring in-depth chemical-sensory analysis of food aroma, as part of the integrated analytical services offered within the RISE research infrastructure.

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