RISE

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Via Orabona 4, 70125 (BA)

ENG: The LeviCell EOS platform enables the enrichment of the target cell population, even from starting cell suspensions with low viability. The platform is based on magnetic levitation technology, in which selected particles are separated through the formation of a density gradient generated by magnets. Magnetic levitation therefore enables cell separation without the use of labels (“label-free” mode), meaning that cells do not undergo alterations to their transcriptional profile caused by binding to the markers used for cell separation in conventional methods, such as FACS and magnetic beads. The system is also based on low-pressure microfluidics (<1 psi), which minimizes induced cellular stress. In the context of single-cell transcriptomic analyses, this technology enables the generation of higher-quality data, as demonstrated by improvements across all key quality metrics, including better median genes/cell, sequencing saturation, fraction of reads/cell, and median UMI counts/cell. The platform uses a single disposable component made of inert and biocompatible plastic, which can accommodate up to four cell samples in parallel. Each assay can process up to 40 million cells, with cell sizes ranging from 1 μm to 350 μm . The target cell population can be readily recovered in 20 minutes, with no risk of cross-contamination between samples. The platform includes both bright-field and fluorescence

microscopy, with two fluorescence channels, allowing cells to be visualized during the separation process. Through the software integrated into the system, cell separation within the cartridge can be monitored in real time, and separation parameters can be adjusted accordingly. The instrument is also equipped with an internal temperature control and regulation system to maintain samples under physiological conditions. The instrument has a compact footprint, allowing it to be placed inside a biological safety cabinet, and does not require calibration or wash cycles between samples.

The platform is equipped with standalone analysis software that allows users to create, visualize, and analyze their samples and is compatible with any PC running Windows 10 or Windows 11. In addition, dedicated software enables run configuration and provides real-time visualization of the entire separation channel for each sample run. During the separation process, the software analyzes the cell population according to density, distribution, and morphology.

Furthermore, the instrument can be equipped with software that enables the acquisition of additional information about samples prior to downstream analysis. Supported by the integrated optical system, including two fluorescence channels, the platform can perform several analytical measurements, such as cell counting, sample quality assessment, debris quantification, and morphological measurements.

Finally, up to four modules can be integrated and connected to a central PC. Each module can be operated and configured independently, thereby increasing overall cell-separation throughput.

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