RISE

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Via Orabona 4, 70125 (BA)

ENG: The iScan system uses its own microarray technology for a wide range of DNA analysis applications, including genome-wide genotyping of humans, plants, livestock, and model organisms. For organisms that have not yet been genetically characterized, customized microarrays can be created using known or putative markers. The iScan system is therefore ideal for rapid, accurate, and high-throughput screening in clinical research, genomics, and agrigenomics applications.

Illumina microarray technology, also known as BeadArray technology, uses silica microspheres housed in precisely etched microwells and coated with hundreds of thousands of copies of a 50-base oligonucleotide probe corresponding to a specific genomic locus. The probe length (50 nucleotides) ensures high hybridization specificity with the target sequence. On the surface of each array, or BeadChip, from hundreds of thousands to millions of genotypes can be analyzed simultaneously for a single individual. When DNA fragments are transferred onto the BeadChip, each probe hybridizes to its corresponding complementary sequence in the DNA sample, leaving the base at the polymorphic site of interest unpaired. A single-base extension performed by a high-fidelity DNA polymerase, incorporating one of four fluorophore-labeled nucleotides, enables detection of the site of interest. When excited by a laser, the fluorophore emits a signal that provides information about the allelic ratio at that

locus. In fact, each assay genotypes a locus using two color channels: one color for each allele. The relative intensity of the two colors indicates whether the genotype at that locus is heterozygous or homozygous.

Genotyping analysis allows the analysis of between 300 and 5 million loci in a single reaction, with specificity provided by hybridization to the 50-mer probe, which ensures specificity at two levels, followed by high-fidelity enzymatic extension.

For CNV analysis, iScan can detect many types and sizes of structural variation in the human genome that affect phenotypes, including duplications, deletions, amplifications, loss of heterozygosity (LOH), and mosaicism. This is possible because dedicated microarrays provide dense coverage of approximately 250 clinically relevant regions, including subtelomeric and pericentromeric regions and the sex chromosomes, which are commonly screened in cytogenetics laboratories. This scalable solution offers substantially higher resolution for detecting smaller regions compared with fluorescence in situ hybridization (FISH) or comparative genomic hybridization (CGH).

The iScan system also enables quantitative measurement of cytosine methylation at individual CpG sites in bisulfite-treated genomic DNA through the use of methylation arrays, providing high resolution for the investigation of epigenetic changes. Following bisulfite treatment, unmethylated cytosine bases are converted to uracil, whereas methylated cytosine bases remain unchanged. The assay interrogates these chemically differentiated loci using two site-specific probes associated with each type of microsphere: one for the methylated locus and the other for the unmethylated locus. The methylation level at the interrogated locus can be determined by calculating the ratio of fluorescent signals from the methylated and unmethylated sites following single-base extension of the probes with a fluorescence-labeled ddNTP.

The iScan system provides analytical flexibility because it uses BeadChips in different formats to analyze 4, 8, 24, or 96 samples in parallel. The system employs high-performance lasers, optics, and detection systems to provide sub-micron resolution and high throughput. Specifically, the instrument is equipped with lasers operating at excitation wavelengths of 532 nm and 658 nm and a scanner with a true spatial resolution of 0.53 microns, featuring dual-channel fluorescence detection, with an analysis time of 20 to 45 minutes per chip. Scanning times are extremely fast without compromising data quality or reproducibility. Samples are prepared manually, but the iScan is compatible with an automated liquid-handling system, which reduces sample preparation time and enables continuous microarray scanning 24 hours a day. The iScan achieves high call rates (>99%), supporting population screening studies and high-resolution CNV analysis, including the accurate detection of single-copy-number variations.

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