
RISE

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Via Edoardo Orabona, 4, 70126,
Bari

ENG: The Nikon Eclipse Ti2 is an inverted fluorescence microscope designed to follow living cells over time with high optical stability and sensitivity. It provides a controlled environment for observing dynamic intracellular events, making it particularly suitable for experiments that require continuous monitoring of calcium signals in the cytosol or mitochondria. When used with ratiometric probes such as Fura-2, the system alternates excitation at the characteristic wavelengths of the dye to generate a fluorescence ratio that reflects changes in intracellular Ca^{2+} concentration. This ratiometric strategy is essential because it compensates for common experimental variability, including uneven dye loading, photo-bleaching, dye leakage, differences in cell thickness, and fluctuations in fluorophore concentration, ensuring that the measured signal truly represents biological changes rather than technical artefacts. The Ti2 supports long time-lapse recordings, multi-position imaging, and automated focus stabilization, allowing researchers to capture calcium transients, oscillatory patterns, and stimulus-evoked responses with both temporal precision and spatial consistency. Its optical configuration and high-numerical-aperture objectives enable quantitative measurements at the level of single cells as well as larger cell populations, maintaining signal fidelity throughout extended acquisitions. The Ti2 is widely used in pharmacological studies, analyses of intracellular signalling pathways, investigations of

mitochondrial function, and cellular toxicity testing. It is a versatile platform for probing calcium dynamics in both excitable and non-excitable cells, offering researchers a reliable window into fast, subtle, and biologically meaningful intracellular events.

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