RISE

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Via Edoardo Orabona, 4, 70125, Bari

ENG: The SOTAX CE 7smart Offline is a flow-through cell dissolution testing system compliant with USP Apparatus 4. It maintains controlled hydrodynamic conditions while automatically collecting samples in tubes or vials, supporting reproducible dissolution and drug-release studies from development through quality control.

The system can operate in open-loop mode, providing sink conditions and automated pH changes, or in closed-loop mode with a defined medium volume. The CP pump delivers a constant flow through the cells; samples are filtered in the flow-through cell head and directed to a C613 or C615 autosampler. An automated cleaning procedure simplifies operation between tests.

The CE 7smart Offline is suitable for active pharmaceutical ingredients, intermediates and finished dosage forms, including oral, topical and parenteral products. It supports biorelevant methods, extended- and delayed-release testing, poorly soluble BCS class II and IV compounds, and in-vitro/in-vivo correlation studies. Methods, results and reports can be managed with WinSOTAXplus software.

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